




P.O. Box 3867
Ithaca, NY 14852 USA
info@mitegen.com

Orders: (607) 266 8877
Tech Support: (607) 592 8812
Fax: (607) 697 0400


MicroTools Kit 2:

20 soft, flexible and durable tools to simplify your sample handling tasks




MicroSievesTM

Sieve, lift, transfer and soak.




MicroSawsTM

Cut protein skins and dense lipid phases




MicroLassosTM

Lasso and pull to detach from surfaces.



MicroTipsTM

Poke, separate, lift and reposition on a MicroMount



MicroSpoonsTM

Push, detach from surfaces, separate from a cluster, lift, transfer. Push down to flatten tip and slide under sample.



MicroChiselsTM

Push, detach from surfaces, separate from a cluster, lift, transfer. Push down to flatten tip and slide under sample.

These tools are thicker and more durable than those in Kit 1, but are less suitable for holding crystals during diffraction experiments.