




P.O. Box 3867  
Ithaca, NY 14852 USA  
[info@mitegen.com](mailto:info@mitegen.com)

Orders: (607) 266 8877  
Tech Support: (607) 592 8812  
Fax: (607) 697 0400

## MicroTools™ Kit 3


20 soft, flexible tools for linear, radial and angular measurements  
Available in both standard and X-ray fluorescent versions



1-5

### Horizontal MicroRuler™


500 µm wide with 25 µm increments  
For linear measurements



6-10

### Vertical MicroRuler™


1000 µm wide with 25 µm increments  
For linear measurements



11-15

### Large MicroProtractor™

1000 µm wide with 50 µm and 15° increments  
For radial and angular measurements



16-20

### Small MicroProtractor™

320 µm wide with 20 µm and 45° increments  
For radial and angular measurements